



### **Laundering and Care**

1. Remove all panels and/ or any metal components prior to cleaning.
2. This brace can be hand washed in cold or lukewarm water with mild soap or detergent and air dried
3. Do not use any substance containing Chlorine

### **Warning**

1. Carefully read this application instruction and warning prior to fitting braces, failure to properly position the brace and fasten each strap will compromise performance and support.
2. This brace will Not prevent or reduce all injuries. Proper rehabilitation and activity modification are essential while wearing the brace.
3. Federal Law restricts the sale of this Device to anyone other than a Licensed Health Care Provider.

### **Caution**

1. Single patient use only
2. Do not over-tighten closures
3. If irritation, pain or discomfort occurs while wearing this brace, discontinue use and immediately consult your medical professional.

EUROPEAN AUTHORIZED REPRESENTATIVE



Wellkang Ltd  
Suite B, 29 Harley Street,  
LONDON, W1G 9QR, U.K.  
[www.CE-Marking.EU](http://www.CE-Marking.EU)



Council Directive  
93/42/EEC of  
14 June 1993  
concerning Medical Device

MANUFACTURED BY

**Dicarre LLC**

[info@dicarre.com](mailto:info@dicarre.com)


**Toll Free: 1-844-DICARRE**

125 Half Mile Road, Suite 200  
Red Bank, NJ 07701

**[www.dicarre.com](http://www.dicarre.com)**

**DICARRE**  
Classique TLSO

DB22  
User Guide  
Application Instruction



1. Always wear a T-shirt or orthopedic body sock under your Dicarre brace.
2. Place the Dicarre orthosis around your waist with side closure attachment on the left side pre-attached to left side of front panel (Fig. 1). Proceed to attach right Velcro placement to right side of front panel (Fig. 2).
3. Align the brace so that the front panel is centered on your abdomen (Fig. 3) and the back panel curve is aligned with the lumbar section of your back (Fig. 4).
4. Detach the elastic Velcro compression straps from the orthosis and pull straps simultaneously from your body (Fig.5).
5. Secure the right and left compression strap to the center of the orthosis (Fig.6).
6. Adjust the thoracic extension to patients height with allen key provided (Fig. 7).
7. Shorten or lengthen shoulder straps as needed to patient.
8. Shoulder straps can be trimmed to fit patients.
9. Finally, ensure that orthosis is securely snug and does not cause any discomfort.

