



Laundering and Care

1. This brace can be hand washed in cold or lukewarm water with mild soap or detergent and air dried
2. Do not use any substance containing Chlorine


Warning

1. Carefully read this application instruction and warning prior to fitting braces, failure to properly position the brace and fasten each strap will compromise performance and support.
2. This brace will Not prevent or reduce all injuries. Proper rehabilitation and activity modification are essential while wearing the brace.
3. Federal Law restricts the sale of this Device to anyone other than a Licensed Health Care Provider.

Caution

1. Single patient use only
2. Do not over-tighten closures
3. If irritation, pain or discomfort occurs while wearing this brace, discontinue use and immediately consult your medical professional.

EUROPEAN AUTHORIZED REPRESENTATIVE



Wellkang Ltd
Suite B, 29 Harley Street,
LONDON, W1G 9QR, U.K.
www.CE-Marking.EU



Council Directive
93/42/EEC of
14 June 1993
concerning Medical Device

MANUFACTURED BY

Dicarre LLC

info@dicarre.com


Toll Free: 1-844-DICARRE

125 Half Mile Road, Suite 200
Red Bank, NJ 07701

www.dicarre.com

DICARRE
Airmesh Pocket™ Sling

DE05
User Guide
Application Instruction



Airmesh Pocket™ Sling

1. Placement

Slide arm into sling as shown.

2. Secure shoulder strap

Wrap shoulder strap around neck and shoulder, position the shoulder pad right on the shoulder, feed the shoulder strap through D-ring and secure.

3. Secure waist strap

Wrap around the waist from the back, feed the waist strap through the D-ring and secure.

4. Fit adjustment

Readjust the strap is necessary ensure the fit is snug and comfortable

