



Laundering and Care

1. Remove all panels and/ or any metal frames or components prior to cleaning.
2. This brace can be hand washed in cold or lukewarm water with mild soap or detergent and air dried
3. Do not use any substance containing Chlorine

Warning

1. Carefully read this application instruction and warning prior to fitting braces, failure to properly position the brace and fasten each strap will compromise performance and support.
2. This support will Not prevent or reduce all injuries. Proper rehabilitation and activity modification are essential while wearing the brace.
3. Federal Law restricts the sale of this Device to anyone other than a Licensed Health Care Provider.
4. Settings of valgus/varus controls should be undertaken by medical professionals only

Caution

1. Single patient use only
2. Do not over-tighten closures
3. If irritation, pain or discomfort occurs while wearing this brace, discontinue use and immediately consult your medical professional.

EUROPEAN AUTHORIZED REPRESENTATIVE



Wellkang Ltd
Suite B, 29 Harley Street,
LONDON, W1G 9QR, U.K.
www.CE-Marking.EU



Council Directive
93/42/EEC of
14 June 1993
concerning Medical Device

MANUFACTURED BY

Dicarre LLC

info@dicarre.com


Toll Free: 1-844-DICARRE

125 Half Mile Road, Suite 200
Red Bank, NJ 07701

www.dicarre.com

DICARRE
Dicarre Challenger Beyond

DK18b
User Guide
Application Instruction



1. Brace is to be worn directly against dry skin, while seated, flex knee to 45 degree and place brace on leg. Fasten the tibia strap as shown in fig 1, position the centre hinge about 1 inch above the knee cap with brace slightly behind midline of the leg as shown in fig 2
2. Starting from the thigh, fasten the straps in proper order as shown in fig 3, 4, 5, 6
3. Fasten the remaining straps.
4. Once you have brace on your leg, walk around for a minute or get the blood stimulated. Sit back down and tighten the straps in proper order as shown in fig 7
5. Hinge stops are provided with brace to control Flexion/Extension as needed for patients as shown in fig 8

