



Laundering and Care

1. This brace can be hand washed in cold or lukewarm water with mild soap or detergent and air dried
2. Do not use any substance containing Chlorine


Warning

1. Carefully read this application instruction and warning prior to fitting braces, failure to properly position the brace and fasten each strap will compromise performance and support.
2. This brace will Not prevent or reduce all injuries. Proper rehabilitation and activity modification are essential while wearing the brace.
3. Federal Law restricts the sale of this Device to anyone other than a Licensed Health Care Provider.

Caution

1. Single patient use only
2. Do not over-tighten closures
3. If irritation, pain or discomfort occurs while wearing this brace, discontinue use and immediately consult your medical professional.

EUROPEAN AUTHORIZED REPRESENTATIVE



Wellkang Ltd
Suite B, 29 Harley Street,
LONDON, W1G 9QR, U.K.
www.CE-Marking.EU



Council Directive
93/42/EEC of
14 June 1993
concerning Medical Device

MANUFACTURED BY

Dicarre LLC

info@dicarre.com

Toll Free-1-844-DICARRE

125 Half Mile Road, Suite 200
Red Bank, NJ 07701

www.dicarre.com

A line drawing of the DW12b Wrist Splint. It shows a rectangular device with a strap and buckle, designed to fit around a wrist. The drawing is simple, using only black outlines on a white background.

DICARRE


**Universal Quick-fit™
Wrist Splint**

DW12b

**User Guide
Application Instruction**

1. Position the Brace

Place the hand into the brace with the thumb inside the thumb hole. The palmer stay should be resting on the palmer section of the hand.



2. Secure Straps

Beginning at the hand and working up the arm, feed each strap through the corresponding D-ring, wrap around hand/wrist.

3. Fit adjustment

Readjust the straps for firm support, the brace should be snug but comfortable.

