



### Laundering and Care

1. This brace can be hand washed in cold or lukewarm water with mild soap or detergent and air dried
2. Do not use any substance containing Chlorine


### Warning

1. Carefully read this application instruction and warning prior to fitting braces, failure to properly position the brace and fasten each strap will compromise performance and support.
2. This brace will Not prevent or reduce all injuries. Proper rehabilitation and activity modification are essential while wearing the brace.
3. Federal Law restricts the sale of this Device to anyone other than a Licensed Health Care Provider.

### Caution

1. Single patient use only
2. Do not over-tighten closures
3. If irritation, pain or discomfort occurs while wearing this brace, discontinue use and immediately consult your medical professional.

EUROPEAN AUTHORIZED REPRESENTATIVE



Wellkang Ltd  
Suite B, 29 Harley Street,  
LONDON, W1G 9QR, U.K.  
www.CE-Marking.EU



Council Directive  
93/42/EEC of  
14 June 1993  
concerning Medical Device

MANUFACTURED BY

**Dicarre LLC**

info@dicarre.com

**Toll Free 1-844-DICARRE**

125 Half Mile Road, Suite 200  
Red Bank, NJ 07701


**www.dicarre.com**



DICARRE  
Sports Ankle Lacer


A line drawing of a right foot wearing a brace. The brace is a wrap-around device with multiple straps and buckles, designed to support the ankle and foot. The drawing is simple, using black outlines on a white background.

DA33  
User Guide  
Application Instruction




### 1. Position the brace

Slip foot into the brace. Ensure the tongue is between the top of the foot and the laces.




### 2. Tighten laces

Pull speed laces tabs to secure.



### 3. Invert foot

Invert foot and wrap the upper lateral strap in a figure 8 around the foot and pull vertically until firm medial tension is achieved. Secure by velcro closure.



### 4. Evert foot

Evert foot and wrap the upper medial strap in a figure around the foot and pull vertically until firm lateral tension is achieved. Secure by velcro closure.



### 5. Wrap and secure calf strap

Wrap upper strap around the calf and secure by velcro closure.